

PRODUCT CATALOGUE CATALOGUE PRODUITS



**Your “One-Stop Shop” for
Healthcare Requirements for**
✓Detection ✓Risk Assessment
✓ Diagnosis ✓Therapy

TRIMPEKS

new health dimensions

Trimpeks Healthcare has been established in 1990 in Istanbul/Turkey, by partners having the awareness of how valuable and important Health is, in a person's life.

Considering that a person should both manage to check himself and/or heal himself, at a daily routine or at the required times; we have built up a complete "Home Healthcare" range.

Considering as well that a person has at the required times, to be cared by the Health Professionals in Hospitals and Clinics; we have also built up a complete "basic" "Hospital Products" range.

In this way, we have become a "**1-Stop-Shop**" for healthcare requirements; which prestigiously serves over 50 countries around the world.

Our criteria for choosing the elements which build up our Product Portfolio is that:
"Each & every Trimpeks Product, should be of quality that can be used with our Family Members"

Since you have this Catalogue in your hand, we welcome you to be a Trimpeks Partner in your country.

Sincerely yours,

Your "TRIMPEKS TEAM"



plusmed

health management



DIABETES CARE.....	4 - 8
TEMPERATURE MEASUREMENT.....	9
BLOOD PRESSURE MEASUREMENT.....	10 - 11
ACCESSORIES.....	12 - 13
STETHOSCOPES.....	14
NEBULIZERS.....	15
OXYGEN CONCENTRATORS.....	16 - 17
PULSEOXIMETERS.....	18
PAP DEVICES.....	19 - 21
SUCTION UNIT.....	22
CPAP MASKS	23
PRESSURE MATTRESS SYSTEMS.....	24 - 25
TENS ELECTRICAL STIMULATOR.....	26
PATIENT CARE.....	27 - 29

pM-100 *fastTest*

Blood Glucose Monitor Device

Recognized for its ease of use and quick response time, the Plusmed *fastTest* Blood Glucose Monitoring System is designed to process accurate results utilizing only a small amount of blood for in vitro diagnostics. The Plusmed *fastTest* Blood Glucose Monitoring System is intended for use in the home and in professional settings to monitor whole blood glucose levels obtained from the fingertip, palm and forearm.




- | | |
|-------------------------|---|
| 1 5 seconds measurement | 7 Big LCD screen |
| 2 0.5 µl blood sample | 8 Hypoglycemic & Hyperglycemic indicator alarms |
| 3 Auto- code system | 9 7-14-30 days average results |
| 4 4 reminder alarms | 10 "Hi" and "Lo" range settings |
| 5 300 memories | 11 Alternative site testing |
| 6 Transfer to PC | |

Order Code: M020751
Product No:8698864011281

pM-100 *fastTest*

Blood Glucose Test Strips

The Plusmed *fastTest* Blood Glucose Test Strips are used with the plusmed *fastTest* Blood Glucose Meter for quantitatively measuring glucose (sugar) in whole blood obtained from the fingertip, palm and forearm. The Plusmed *fastTest* Test Strips are for testing outside the body (in vitro diagnostic use).



- 1 For testing glucose in capillary blood
- 2 Tiny blood sample size (0.5 µl)

Order Code: M020752
Product No:5060186753621

fasTest DailyCheck

Blood Glucose Monitoring System


Recognized for its ease of use and quick response time, the Plusmed fasTest DailyCheck Blood Glucose Monitoring System is designed to process accurate results utilizing only a small amount of blood for in vitro diagnostics. The Plusmed fasTest DailyCheck Blood Glucose Monitoring System is intended for use in the home and in professional settings such as hospitals and clinics, to monitor whole blood glucose levels obtained from alternate sites in different positions such as sitting on a chair or at bedside. The display with contrast, large icons and large numbers enables people to read the test results easily.



Order Code: M020775

Product No:8698864016071

- 1 5 seconds measurement
- 2 1 µl blood sample
- 3 No coding required
- 4 Alarm for testing schedule
- 5 360 test memories with date&time
- 6 7, 14 and 28 day averaging
- 7 mmol/L or mg/dL (switchable)



Order Code: M020776

Product No:8698864016088

pM-900 fasTest DiagnoPRO

Blood Glucose Monitoring System

Recognized for its ease of use and quick response time, Plusmed fasTest DiagnoPro Multi-Monitoring System is designed to process accurate results utilizing only a small amount of blood for in vitro diagnostics. Plusmed fasTest DiagnoPro Multi-Monitoring System is for quantitatively measuring Blood Glucose (sugar), Blood -Ketone and Blood Uric Acid in fresh capillary whole blood, neonatal blood, arterial and venous blood samples. It is only for use outside the body (For in vitro diagnostic use) and intended for self-testings of patients at home and professional use in hospitals and clinics to monitor blood glucose levels.



Order Code: M020721

Product No:8698864014541

- 1 Multi Function (Blood Glucose, Blood - Ketone* and Blood Uric Acid*)
- 2 GOD 0.7 µl Blood Sample
- 3 Average glucose in 7, 14, 30, 60 and 90-days
- 4 Results in 5 Seconds
- 5 Memory Up to 500 Tests
- 6 Strip Ejector
- 7 ISO 15197:2013 / EN ISO 15197:2015 Certified
- 8 PC Download via data transfer cable
- 9 Suitable for Neonatal use
- 10 6 months test strip expiry after vial's first opening
- 11 Measurement range indicator



Order Code: M020722

Product No:8698864014572

fastest Hblyzer

HB Hemoglobin Testing System

Fasttest Hblyzer is intended to be used for quantitative determination of Hemoglobin (Hb) and calculated Hematocrit (Hct) in capillary and venous human whole blood. This system is intended for outside of the body (in vitro diagnostic use) by physical examination and by healthcare professional in a clinical situation such as anaemia.

- 1 Results in 15 seconds
- 2 10 ul small blood sample size
- 3 Memory stores 1000 reading with date and time
- 4 Switchable three measuring unit
- 5 Deluxe carry case
- 6 Large LCD
- 7 Auto power off
- 8 Uploading Data to PC via Plusmed Fasttest Hblyzer Software



Order Code: M020783

Product No: 8698864014954

fastest Hblyzer

HB Hemoglobin Test Strips

fastTest Hblyzer hemoglobin strip is used with fastTest Hblyzer hemoglobin meter as a Hemoglobin Analyzing System to quantitatively measure hemoglobin with capillary and venous whole blood. The Hemoglobin Analyzing System is intended for in vitro diagnostic use (for Medical Detection Mechanism) only (outside of the body).

- 1 For testing hemoglobin with capillary and venous whole blood
- 2 Tiny blood sample volume (10 μ L)
- 3 50 test strips



Order Code: M020781

Product No: 8698864014961

GLYCOHEMOGLOBIN ANALYZER

Unique Handheld HbA1c. Analyzer for Better Diabetes Care

- 1 Portable Testing of HbA1c
- 2 Accurate Result with CV<3%
- 3 Easy and Fast Operation
- 4 Room Temperature Storage



P3 (CGM)

Real-Time Continuous Glucose Monitoring System

Accurate and Affordable Continuous Glucose Monitoring System

- 1 14 days life
- 2 MARD=10%
- 3 Comfortable
- 4 Real Time
- 5 Algorithm
- 6 Cloud Portal



fasTtest

Blood Lancets

Plusmed fasTTest blood lancets are used with lancing device for taking blood sample. Plusmed fasTTest blood lancets are disposable and for single patient use.



- 1 Painless application due to thinness
- 2 Needle gauge: 30G
- 3 More reliable due to gamma sterilization method

Order Code: M021962
Product No:8698864013575

FASTFINE

Insulin Pen Needles

Pen needles are used in conjunction with injection pens to deliver injectable medications into the body. Insulin pen needles come in a variety of lengths and widths to suit all body types. From 8 mm to 4 mm, the needle length you choose is likely to be dictated by your size, children being likely to benefit from the shorter 4 mm size. Each pen needle is embedded in a plastic hub that can attach to the pen.



4MM-32G-100
Order Code: M026605
Product No:6945630106408

5MM-32G-100
Order Code: M026609
Product No:6945630106613

6MM-32G-100
Order Code: M026606
Product No:6945630106385

8MM-31G-100
Order Code: M026607
Product No:6945630106392

- 1 Sterile
- 2 Fits all insulin pens
- 3 Non toxic, non pyrogenic
- 4 Silicon coated
- 5 Available in 4 different sizes
 - 32 G (0.23 X 4 mm)
 - 32 G (0.23 X 5 mm)
 - 31 G (0.25 X 6 mm)
 - 31 G (0.25 X 8 mm)

pM-101

Pen Type - Digital Thermometer

This thermometer can measure the temperature in various locations on the body which maintain a fairly stable temperature (oral, axillary, rectally). Display readings to great precision (0.1 °C or 0.2 °F).

- 1 ± 0.1 °C Accuracy
- 2 Fast response
- 3 Memory function
- 4 Auto shut-off



Order Code: M029421

Product No:8698864013049

UltraScan II

Infrared Thermometer
Ear & Forehead Scan

- 1 Measured on forehead or ear
- 2 Fast reading
- 3 10 reading memory
- 4 3-color Backlight
- 5 Result in 1 second
- 6 Automatic power-off



Order Code: M020438

Product No:6970392214481

VOICE

Arm-Type Fully Automatic TALKING
Digital Blood Pressure Monitor

- 1 Accurate Measurement
- 2 Quiet Inflation Method
- 3 Automatic Power Off
- 4 Irregular Heartbeat Indicator
- 5 Last 3 Tests Average
- 6 Voice
- 7 Static Pressure Measurement
- 8 2x150 Memories in Two Groups with Date and Time
- 9 Cuff Circumference Approx.135 (W) x 485(L) mm

Order Code: M028742
Product No:6970392216799



CONNECT

Arm-Type Fully Automatic
Digital Blood Pressure Monitor

- 2 Accurate Measurement
- 3 Quiet Inflation Method
- 3 Automatic Power Off
- 4 Irregular Heartbeat Indicator
- 5 Last 3 Tests Average
- 6 Static Pressure Measurement
- 7 Bluetooth connection
- 8 2x150 Memories in Two Groups with Date and Time
- 9 Cuff Circumference Approx.135 (W) x 485(L) mm

Order Code: M028743
Product No:6970392213408



pM-8197W

Wrist-Type Fully Automatic Digital Blood Pressure Monitor

- 1 Average of Last 3 Measurements
- 2 Accurate Measurement
- 3 Irregular Heartbeat Detection
- 4 Quick Result
- 5 Quiet Inflation Method
- 6 Automatic Power-off
- 7 Wrist Circumference: 13,5-21,5 cm (5,3-8,5")

Order Code: M028740
Product No:8698864016484



PERFECT

Professional Aneroid Sphygmomanometer

- 1 Chrome plated metal case
- 2 48 mm diameter of dial
- 3 Stainless steel handle with chrome plated
- 4 0-300 mmHg scale
- 5 ± 3 mmHg measurement accuracy
- 6 Green Cotton Adult Hook Cuff
- 7 Double tube system
- 8 Storage case

Order Code: M021629

Product No:6970384252132



Palm Sphyg I

Aneroid Sphygmomanometer

- 1 ABS plastic upper case
- 2 60.2 mm diameter of dial
- 3 Right or left hand usage
- 4 Large zinc-alloy handle with black powder-coating
- 5 Deflation valve: "Push to button" for rapid deflation, "screwing" for precise deflation
- 6 Nylon D-Ring Adult Cuff
- 7 Silicone shockproof ring
- 8 Storage case



Order Code: M021630

Product No:6970384252187

pM-A01S

Aneroid Sphygmomanometer

Blood pressure measurements with the PM-A01S are done using the cuff / stethoscope listening method. It provides more reliable measuring.

- 1 300 mmHg non-stop pin manometer
- 2 D-ring nylon cuff
- 3 Heavy duty metal gauge
- 4 Lightweight exterior single head stethoscope
- 5 Chrome plated brass air relief valve



Order Code: M029330

Product No:8698864013230

METAL HOOK CUFF

- 1 Compatible with Aneroid Sphygmomanometers
- 2 Double tube
- 3 Latex bladder
- 4 Green
- 5 Cuff dimensions: 49 x 14 cm
Bladder dimensions: 22 x 12 cm
Arm circumferences: 25 – 40 cm



MANOMETER

Compatible with Aneroid Sphygmomanometers



BULB

Compatible with Aneroid Sphygmomanometers



VELCRO CUFFS



Specification	Size (mm)	Arm Range (mm)
New born	185*55	---
Infant	255*75	130 ~ 190
Child	345*110	184 ~ 267
Adult	540*145	254 ~ 406
Large adult (Obese)	620*175	333 ~ 510
Thigh	780*220	406 ~ 660

D-RING CUFFS



Specification	Size (mm)
Adult	500*140
Large Adult (Obese)	600*170

pM-ST71

Stethoscope - Single Head

- | | |
|---------------------------------------|---|
| 1 Comfortable fit | 6 Aluminum alloy |
| 2 High Quality | 7 Aluminum alloy binaural with chromium plate |
| 3 Better healthcare | 8 Chest piece diameter: 46 mm |
| 4 Good acoustic response | 9 Diaphragm diameter: 42.8 mm |
| 5 Ear tips with soft cone-shaped bell | |



Order Code: M020180
Product No: 8698864015524

pM-ST72

Stethoscope - Dual Head

- | | |
|-------------------------------|---|
| 1 Comfortable fit | 6 Ear tips with soft cone-shaped bell |
| 2 High Quality | 7 Opposite direction, two feature function |
| 3 Better healthcare | 8 Aluminum alloy |
| 4 Good acoustic response | 9 Aluminum alloy binaural with chromium plate |
| 5 Bell with non-chilling ring | 10 Chest piece diameter: 46 mm |
| | 11 Diaphragm diameter: 42.8 mm |



Order Code: M020186
Product No: 8698864015531

CardioSTETH I

Stainless Steel Cardiology Stethoscope

- 1 Silver dual head chestpiece for cardiology diagnosis
- 2 Stainless steel chestpiece, excellent sound conduction
- 3 The membrane is sensitive for sound collecting
- 4 A" shaped thick tube with excellent resistance to ambient interference ensuring best sound conduction
- 5 Comes complete with accessories



Order Code: M020190
Product No: 6970384252194

COMPACT NEB

Piston Compressor Nebulizer

- 1 Ideal for All Ages
- 2 Effective medication nebulization
- 3 Handle for Easy Carry
- 4 One-Button Operation
- 5 Light Weight
- 6 Easy to Clean
- 7 Latex Free



Order Code: M029560
Product No: 8698864016798

RespiAid Mesh-01

Mesh Nebulizer

Plusmed RespiAid Mesh-01 is intended for the treatment of respiratory diseases such as asthma, allergies, etc. Thanks to the self-excited oscillation circuit, the ceramic sheets contained in nebulizing head together with the metal mesh are driven successively to oscillate at a high frequency, which makes the liquid flowing through the micromesh of the metal mesh and being ejected and nebulized to form an inhalant aerosol.

- 1 Replaceable medicine cup
- 2 Quiet design
- 3 Dual power supply (2 x 1.5 V AA batteries or Micro USB 5 V)
- 4 Particle size $\leq 5\mu$
- 5 Nebulization rate ≥ 0.20 ml/min
- 6 Capacity of medicine cup ≤ 8 ml
- 7 Accessories: Medicine Cup | Mouthpiece | Adult Mask | Child Mask | USB Cable | Storage Bag | Battery



Order Code: M029557
Product No: 8698864015531

RespiAid P2

Portable Oxygen Concentrator

Plusmed RespiAid P2 Portable Oxygen Concentrator is used on a prescriptive basis by patients requiring supplemental oxygen. It supplies a high concentration of oxygen to the patients.

It supplies a high concentration of oxygen and is used with a nasal cannula to channel oxygen from the concentrator to the patient. The Plusmed RespiAid P2 is small, portable and may be used in home, institution and various mobile environments. Plusmed RespiAid P2 Portable Oxygen Concentrator conforms to all applicable FAA acceptance criteria for Portable Oxygen Concentrator carriage and use on board aircraft.

- 1 Lightweight 1.98kg
- 2 Mini Size: 221*85*160mm, Very small size, easy to carry
- 3 5 Setting Flows: From 1-5 setting flows, can meet different oxygen demand
- 4 Pulse Oxygen Supply: Plusmed RespiAid P2 is equivalent to 5 liters
- 5 High Efficiency Molecular Sieve: Efficient adsorption capacity, can purify more pure
- 6 Oxygen Concentration: 90% -3% / +6%
- 7 Powerful Compressor: Powerful compressor can provide stable output
- 8 Triple Filtration: Effectively remove air pollen, large particles, dust and other impurities
- 9 Battery: 8 cell standard Battery Up to 4 hours



Order Code: M029557
Product No: 8698864016378

pM-KN01

Stationary Oxygen Concentrator

pM-KN01 Oxygen Concentrator, adopt pressure swing adsorption principle, which can separate oxygen, nitrogen and other gas from the air. As soon as power is connected, the oxygen that meets medical use standards can be separated from air constantly. Planned life of the whole generator reaches up to 20,000 hours. There is no influence on indoor oxygen percent during the generator operating.

- ① Oxygen Concentration: 93%±3%Oxygen
- ② Flow: 0-5 L/min
- ③ Outlet Pressure: 0.04-0.07 Mpa
- ④ LCD Display
- ⑤ Alarms: Power failure alarm, Low purity alarm, High pressure alarm, Low pressure alarm
- ⑥ Fire Safe Valve

Order Code: M023022
Product No:8698864012950



pM-KN10

Stationary Oxygen Concentrator

- ① Flow rate: 0-10L/min
- ② Oxygen concentration: 1-9L/min: 93% (± 3%); 10L/min: 90% (± 2%)
- ③ Fire Safe Valve
- ④ Outlet pressure: 0.04-0.07Mpa
- ⑤ Noise level: ≤50dB
- ⑥ Power: AC220V (±10%), 50/60Hz (±1Hz); AC110V (±10%), 50/60Hz (±1Hz)
- ⑦ Power consumption: ≤550W
- ⑧ LCD display
- ⑨ Alarms: Power failure alarm, Low purity alarm, High pressure alarm, Low pressure alarm
- ⑩ Net weight: 28kg
- ⑪ Size: 365×375×600mm

Order Code: M023024
Product No:8698864015500



OxyCheck

Pulseoximeter

OxyCheck Pulse Oximeter is a non-invasive device intended for the spot-check of oxygen saturation of arterial hemoglobin (SpO₂) and the pulse rate of adult and pediatric patients in home and hospital environments (including clinical use in internist/surgery, anesthesia, intensive care etc.).

- 1 Big display
- 2 Measures blood oxygen saturation level and pulse rate
- 3 SpO₂ & PR waveforms
- 4 Perfusion Index (PI)
- 5 Easy to Operate
- 6 High Accuracy
- 7 Low battery indicator
- 8 Automatic power off within 5 seconds
- 9 Screen displayed in 4 directions



Order Code: M061607

Product No:8698864015562

OxyCheck II

Pulseoximeter

Plusmed OXYCHECK II is a kind of innovated medical detection device with non-invasive and continuous features for artery SPO₂ and PR detection. Plusmed OXYCHECK II is portable and easy to measure the SPO₂ and PR value quickly and precisely. Plusmed OXYCHECK II measures human blood oxygen saturation and heart rate through the finger Pulse Oximeter.

Plusmed OXYCHECK II is suitable for family, clinic, oxygen bar, sports health (use before and after exercise. It is not recommended for use during exercise), community health and other ranges. Plusmed OXYCHECK II is not suitable for monitoring the patient's prolonged use.

- 1 Big display
- 2 Measures blood oxygen saturation level and pulse rate
- 3 SpO₂ & PR waveforms
- 4 Perfusion Index (PI)
- 5 Easy to Operate
- 6 High Accuracy
- 7 Low battery indicator
- 8 Automatic power off within 8 seconds
- 9 Screen displayed in 4 directions



Order Code: M061608

Product No:8698864016156

RespiAid

CPAP SYSTEM

RespiAid CPAP is equipped with required alarm and warning systems to offer simplicity, safety and comfort to the user to provide continuous treatment.



- | | |
|--|---|
| <ul style="list-style-type: none"> ① Modes: CPAP ② Pressure Range: 4 – 20 hpa ③ Ramp time: 0-60 minutes ④ LCD Display ⑤ Powerful Analysis Software ⑥ SD card ⑦ Auto Altitude Compensation ⑧ Auto Leak Compensation ⑨ Auto on/off ⑩ High performance motor ⑪ Humidification: 0-5 | <ul style="list-style-type: none"> ⑫ Humidifier: <ul style="list-style-type: none"> Easy clean Easy to fill water chamber Metal Material Fast Heating Good humidification ⑬ Noise Level: <30 dB ⑭ Language: Multiple ⑮ Weight: <ul style="list-style-type: none"> Without humidifier : <1 kg Humidifier included : <1.6 kg ⑯ Dimensions: <ul style="list-style-type: none"> 140 mm × 156mm × 94 mm (without Humidifier) 262 mm × 156mm × 94mm (with Humidifier) |
|--|---|

Order Code: M022507
Product No: 8698864015647

RespiAid

AutoCPAP SYSTEM

Designed to provide the user with a comfortable sleep, RespiAid Auto CPAP offers all kinds of details for tracking the treatment parameters, and as an extremely technological product, it is equipped with alarm and warning systems and provides simplicity, safety and comfort to the user.



- | | |
|---|---|
| <ul style="list-style-type: none"> ① Modes: CPAP, APAP ② Pressure Range: 4 – 20 hpa ③ Ramp time: 0-60 minutes ④ LCD Display ⑤ Powerful Analysis Software ⑥ SD card ⑦ EP Flex ⑧ Auto Altitude Compensation ⑨ Auto Leak Compensation ⑩ Auto on/off ⑪ High performance motor ⑫ Humidification: 0-5 | <ul style="list-style-type: none"> ⑬ Humidifier: <ul style="list-style-type: none"> Easy clean Easy to fill water chamber Metal Material Fast Heating Good humidification ⑭ Noise Level: <30 dB ⑮ Language: Multiple ⑯ Weight: <ul style="list-style-type: none"> Without humidifier : <1 kg Humidifier included : <1.6 kg ⑰ Dimensions: <ul style="list-style-type: none"> 140 mm × 156mm × 94 mm (without Humidifier) 262 mm × 156mm × 94mm (with Humidifier) |
|---|---|

Order Code: M022508
Product No: 8698864015654

RespiAid

BPAP-30T – BPAP S SYSTEM

Plusmed RespiAid BPAP-30T “S BPAP” can be comfortably used during central apnea besides supporting spontaneous breathing. While providing respiration with minimum effort thanks to the improved breathing compatibility technology without tiring the user. Furthermore, it supports the healthcare professionals with detailed analysis and an integrated reporting program.



Order Code: M022509

Product No: 8698864015661

- | | |
|-------------------------------------|--|
| 1 Modes: CPAP, APAP, BPAP-S, BPAP T | 13 Humidifier:
Easy clean |
| 2 Pressure Range: 4 – 30 hpa | Easy to fill water chamber |
| 3 Ramp time: 0-60 minutes | Metal Material |
| 4 LCD Display | Fast Heating |
| 5 Powerful Analysis Software | Good humidification |
| 6 SD card | 14 Adjusting the ventilation parameters in real time |
| 7 EP Flex | 15 Noise Level: <30 dB |
| 8 Auto Altitude Compensation | 16 Language: Multiple |
| 9 Auto Leak Compensation | 17 Weight: |
| 10 Auto on/off | Without humidifier : <1 kg |
| 11 High performance motor | Humidifier included : <1.6 kg |
| 12 Humidification: 0-5 | 18 Dimensions: |
| | 140 mm × 156mm × 94 mm (without Humidifier) |
| | 262 mm × 156mm × 94mm (with Humidifier) |

RespiAid

BPAP-30T Pro – BPAP ST SYSTEM

Plusmed RespiAid BPAP-30T Pro “ST BPAP” can be comfortably used during central apnea besides supporting spontaneous breathing. While providing respiration with minimum effort thanks to the improved breathing compatibility technology without tiring the user. Furthermore, it supports the healthcare professionals with detailed analysis and an integrated reporting program.



Order Code: M022510

Product No: 8698864015678

- | | |
|--|--|
| <ul style="list-style-type: none"> ❶ Modes: CPAP, APAP, BPAP-S, BPAP-T, BPAP-ST ❷ Pressure Range: 4 – 30 hpa ❸ Ramp time: 0-60 minutes ❹ LCD Display ❺ Powerful Analysis Software ❻ SD card ❼ EP Flex ❽ Auto Altitude Compensation ❾ Auto Leak Compensation ❿ Auto on/off ⓫ High performance motor ⓬ Humidification: 0-5 | <ul style="list-style-type: none"> Ⓜ Humidifier: <ul style="list-style-type: none"> Easy clean Easy to fill water chamber Metal Material Fast Heating Good humidification Ⓨ Adjusting the ventilation parameters in real time Ⓩ Noise Level: <30 dB ⓐ Language: Multiple ⓑ Weight: <ul style="list-style-type: none"> Without humidifier : <1 kg Humidifier included : <1.6 kg ⓓ Dimensions: <ul style="list-style-type: none"> 140 mm × 156mm × 94 mm (without Humidifier) 262 mm × 156mm × 94mm (with Humidifier) |
|--|--|

pM-7E A

Electrical Phlegm Suction Unit

Portable absorb phlegm apparatus designed based on developing orientation of similar products at home and abroad is a new generation of oil free lubrication suction device, which is suitable for use by the patient who has difficulty in phlegm removal due to illness, coma and operation, as well as for aspirating such liquid as pus and blood during the clinical practice. It is the commonly applied medical device for use in the emergency room, operation room, and for nursing in sickroom and home health care.



- | | |
|-----------------------------|--------------------------|
| 1 High negative pressure | 4 Embedded liquid holder |
| 2 Large flux | 5 Low noise operation |
| 3 Oil free lubrication pump | |

Order Code: M024901
Product No:8698864001015

pM-AF03

Hydrophobic Suction Filter

Size: H:54mm; w. 64mm

Housing:ABS

Filter: Hydrophobic PTFE Media Disc

BFE: 99.999%

Pore Size:0.2um

Weight:18.2g

Connector:8mm

Stepped Hose Barbs

Max Operating Pressure : 20psi

Effective Filtration Area: 24.6cm²

Max Operating Temp:60°C



pM-AF05

Hydrophobic Suction Filter

Size: H:54mm;w:58mm

Housing:PPFilter:

Air filtration membrane

BFE:99.999%

Pore Size:0.3um

Weight:19.8g

Connector:6mm-8mmStepped Hose Barbs

Max Operating Pressure:20psi

Effective Filtration Area: 15.5cm²

Max Operating Temp: 60°C



pM-FM

Full Face CPAP Mask

pM-FM Face Mask is comfortable flexible and quiet as it is designed to minimize the contact with the skin, which provides a comfortable feeling to the user during the treatment.

- 1 Anti-asphyxia valve
- 2 Adjustable silicone forehead support
- 3 Two sided ventilation
- 4 Oxygen entry port
- 5 Latex free
- 6 Adjustable Clips
- 7 Revolver head elbow



MEDIUM

Order Code: M027181

Product No:8698864011687

LARGE

Order Code: M027183

Product No:8698864011687

pM-NM

Nasal CPAP Mask

pM-NM Nasal Mask is designed for a comfortable, quiet and flexible use. Due to the minimized contact with the skin, the user has a comfortable feeling during the treatment.

- 1 Adjustable silicone forehead support
- 2 Ventilation
- 3 Oxygen entry port
- 4 Does not Contain Latex
- 5 Unbreakable clips
- 6 Revolver head elbow



MEDIUM

Order Code: M027180

Product No:8698864011670

LARGE

Order Code: M027182

Product No:8698864011670

pM-AD01

Tubular Pressure Mattress System

This product is intended to help and reduce the incidence of pressure ulcers while optimizing patient comfort. For emergent CPR, please pull the CPR yellow strip outside quickly to reduce the air fast.

- 1 Air ventilation function
- 2 Tubes on A-B lines provide alternating inflating
- 3 CPR - Fast Air Release Valve
- 4 Made of heavy duty, premium PVC
- 5 Easy to use
- 6 Adjustable
- 7 Low noise, low vibration, low temperature
- 8 Easy fixing on bed frame with hooks

Order Code: M022141

Product No:8698864007123



pM-2012 PLUS

Bubble Pressure Mattress System

This product is intended to help and reduce the incidence of pressure ulcers while optimizing patient comfort. It also provides individual home care setting and long term care of whom suffering from pressure ulcer and pain management as prescribed by a physician.

- 1 Easy to use
- 2 Adjustable
- 3 Low noise
- 4 Low vibration, low temperature

Order Code: M022180

Product No:8698864013322



pM-TE01

TENS + EMS Electrical Stimulator

Plusmed TENS+EMS Dual Therapy System is a dual channel digital TENS/ NMES electrical stimulator for active treatment application. Transcutaneous Electrical Nerve Stimulator (TENS) can be used for symptomatic relief of chronic intractable pain. Neuromuscular electrical stimulation (NMES) can be used for relax muscle spasms, prevent or retard disuse atrophy, increase local blood circulation, re- educate muscles or maintain / increase the range of motion.

- 1 Digital technology
- 2 2 channels
- 3 8 different programs
- 4 Low battery indicator
- 5 12-step-timer
- 6 Adjustable intensity and treatment time

Order Code: M026501

Product No:8698864007703



pM-CH320

Reuseable Hot / Cold Pack

The use of heat and cold is a standard treatment for a variety of injuries. The type and age of the injury or condition indicates whether hot or cold is appropriate. The use of hot and cold packs is one of the simplest methods of applying heat or cold to the affected area. Cold pack decrease the bruising associated with sprains, as well as the swelling and help to relieve the pain associated with these conditions by numbing the area. Hot pack relaxes muscles and improves blood circulation and produces a soothing effect on the body, further aiding in the healing process.



- 1 Comfortable and convenient application
- 2 Provides effective relief from muscular pain
- 3 Helps reduce swelling
- 4 Great for treating a wide range of tissue injuries
- 5 Reusable pack is perfect for continuous treatment at home

Order Code: M024857
Product No:8698864007758

pM-CS01

Local Anesthetic Freezing Spray

pM-CS01 Ethyl chloride is manufactured to make local anesthesia for skin by providing excessive cooling. Further, it can be used for the injuries of the soft tissue, edema and similar illnesses that occur because of sport activities and for small operations like wart removal.



- 1 Anesthesia prior to injection
- 2 Relief in sport injuries
- 3 Eliminates edema

Order Code: M020601
Product No:8698864011267

PULL UPS ADULT PANTS



- 1 360° waistband-close
- 2 Leakproof barriers
- 3 Super absorbent inner surface
- 4 Triple textile surface
- 5 Sizes: Large, Medium

ADULT DIAPER



- 1 Breathable textile surface
- 2 The inner surface absorbs and traps the liquid rapidly
- 3 Leakproof barriers
- 4 Elastic waistband
- 5 Sizes: Large, Medium, Small

DISPOSABLE UNDERPADS

- 1 Skin-friendly soft texture
- 2 Leakage protection edges
- 3 Long lasting backsheet
- 4 Durable outer surface
- 5 60x90 cm size



HYGIC PERINEAL CLEANING TOWEL



- 1 Alcohol-free
- 2 Latex-free
- 3 Lanolin-free
- 4 Anti-fungal
- 5 Anti-allergy
- 6 Anti-bacterial

EXAMINATION TABLE COVER



- Professional subsidiary products eases the medical examination and bedridden patients body cleaning

COLOSTOMY - ILEOSTOMY - UROSTOMY SYSTEM



- 1 One-Piece Colostomy Pouch (Closed or Drainable)
- 2 One-Piece Urostomy Pouch (With Tap)
- 3 Two-Piece Colostomy Pouch (Closed or Drainable)
- 4 Baseplate for Two-Piece Colostomy System
- 5 Baseplate for Two-Piece Urostomy System





vivocare
health enjoyment



Index



THERMOMETER	31
BLOOD PRESSURE MONITOR.....	31
NEBULIZER.....	32
BABY CARE.....	33

TRULY

Digital Thermometer - INFRARED NON-CONTACT

It measures the ambient and body temperature without contact with colored warning symbols.

- 1 0,5 second measurement
- 2 LCD display with backlight
- 3 Fever alarm
- 4 32 memories
- 5 °C or °F switchable
- 6 Automatic shut-off
- 7 Color Temperature Warner

Order Code: M029631

Product No: 8698864013742



GUARDY

Blood Pressure Monitor - AUTOMATIC UPPER ARM

„Guardy“ can be used by 4 users and is monitoring the risk of Hypertension according to WHO.

- 1 Irregular Heartbeat Detector (IHB)
- 2 WHO risk category indicator
- 3 120 memories with date and time for 4 users
- 4 Average of the last 3 measurements
- 5 Auto power off

Order Code: M028721

Product No: 8698864007697



STEAMY

Nebulizer - COMPRESSOR TYPE

STEAMY is a user friendly nebulizer with a powerful compressor. It is designed to deliver the prescribed medication solution to treat patient respiratory disorders, such as asthma, allergies and bronchitis. The nebulizer converts the medication solution into an aerosolized mist which is inhaled by the patient through the mouthpiece or masks.

- 1 Average Nebulization Rate > 0.32ml/min
- 2 Mode of Operation: 24 hours continuous
- 3 Flow Rate: 5-8 LPM
- 4 MMAD: 3 μ m



Order Code: M029724
Product No: 8698864013339

LUCKY

Nebulizer - COMPRESSOR TYPE

LUCKY has an adorable, cute, child-friendly design

- 1 Effective medication delivery.
- 2 Average Nebulization Rate > 0.32ml/min
- 3 Mode of Operation: 24 hours continuous
- 4 Flow Rate: 5-8 LPM
- 5 MMAD: 3 μ m



Order Code: M029723
Product No: 8698864012127

FONDLY

ELECTRONIC BABY SCALE

- 1 4 high precision strain gauge sensor
- 2 Auto on and auto off
- 3 Tare function
- 4 Auto lock/hold function
- 5 Low battery and overload indication
- 6 Unit conversion switch
- 7 With Height Measurement Tape



Order Code: M024324

Product No: 8698864015753

SMILKY

Breast Pump – ELECTRONIC

This electric breast pump is controlled by microcomputer chip. It has 2-phase suction which mimic sucking actions of infant. More milk for the baby, more comfortable for the Mom. The stimulating speed is fixed and stimulating suction divided into 9 adjustable levels. The expressing speed and suction are divided into 9 adjustable levels. It can mimic baby's sucking rhythm as well as sucking quantity. It can save or run the record of sucking speed and suction setting beforehand.

- 1 Adjustable 9 speeds
- 2 Silicone cushion
- 3 BPA free
- 4 LCD Screen
- 5 Memory function: Save the setting of last work under standby
- 6 Portable: Adapter & battery operated with lanyard and hook design
- 7 Quiet design



Order Code: M027501

Product No: 8698864007666



F. BOSCH
MEDIZINTECHNIK



Index



BLOOD PRESSURE MONITORS	35
NEBULIZER.....	36
FETAL DOPPLER.....	37
PAIN THERAPY	38
PHLEGM SUCTION MACHINE	39

PRAKTISCH

Talking Automatic Digital Blood Pressure Monitor

- 1 Big LCD screen
- 2 120 Memories
- 3 2 languages chooseable
- 4 Date and time display
- 5 Motion detection
- 6 Hypertension risk indicator
- 7 Irregular Heartbeat Detector (IHB)
- 8 Average of last 3 readings
- 9 Low battery indicator
- 10 Auto power off



Order Code: M028723
Product No: 8698864015494

INTELLIGENT

Arm-Type Fully Automatic Digital Blood Pressure Monitor

- 1 Average of Last 3 Measurements
- 2 Irregular Heartbeat Detection
- 3 Accurate Measurement
- 4 Quick Result
- 5 Quiet Inflation Method
- 6 Automatic Power-off

Order Code: M028738
Product No: 6970392213477



GENTLE

Wrist-Type Fully Automatic Digital Blood Pressure Monitor

- 1 Average of Last 3 Measurements
- 2 Accurate Measurement
- 3 Irregular Heartbeat Detection
- 4 Quick Result
- 5 Quiet Inflation Method
- 6 Automatic Power-off
- 7 Wrist Circumference: 13,5-21,5 cm (5,3-8,5")

Order Code: M028739
Product No: 6970392214368



MOBINEB

Mesh Nebulizer

- 1 Dual Power supply with choice of Rechargeable Battery or USB Cable
- 2 Two Level Nebulization Rate Adjustment
- 3 Tiny particle size < 5 μm
- 4 Big medication cup
- 5 Low noise
- 6 Light-Weight
- 7 Easy to Use

Order Code: M029559

Product No: 8698864016453



BABYGERÄUSCH

Fetal Doppler

- 1 Pocket type
- 2 Accurate Detection
- 3 Rechargeable Battery
- 4 Lightweight
- 5 LCD Display
- 6 Backlight
- 7 3.0 MHz Ultrasound Frequency
- 8 Battery Indicator
- 9 3 FHR modes: Manual Counting, Real-Time & Average
- 10 Adjustable Volume Level
- 11 Audible Alarm

Order Code: M067713

Product No: 8698864016453



FB2405

TENS + EMS Dual Therapy Electrical Stimulator

- 1 8 different therapy programs
- 2 2 channels
- 3 Low battery indicator
- 4 Adjustable intensity and treatment time
- 5 12-step-timer
- 6 Auto shut-off

Order Code: M029001

Product No:4260157890213




PROVAC

Portable Phlegm Suction Machine

Portable absorb phlegm apparatus designed based on developing orientation of similar products at home and abroad is a new generation of oil free lubrication suction device, which is suitable for use by the patient who has difficulty in phlegm removal due to illness, coma and operation, as well as for aspirating such liquid as pus and blood during the clinical practice. It is the commonly applied medical device for use in the emergency room, operation room, and for nursing in sickroom and home health care.

- 1 With a low flow, medium vacuum and fashionable stylish
- 2 Adopting oil-free piston pump
- 3 Adjustable pressure button to select any suitable pressure level
- 4 Overflow protection can safely and efficiently prevent the liquid from entering into the pump

Order Code: M024906
Product No: 8698864015548





promedic
health treatment



Index



TOURNIQUET SYSTEM.....	41
PATIENT MONITORS.....	42
HOLTER SYSTEM.....	43
HANDHELD PULSEOXIMETER.....	43
ECG BELT.....	43
ELECTROSURGICAL UNITS.....	44 - 45
FETAL MONITORS.....	46
LARYNGOSCOPE SETS.....	47
DIAGNOSTIC SETS.....	47 - 48

DTS-3000

DIGITAL TOURNIQUET SYSTEM

- 1 Separate pressure and operation for 2 channels
- 2 Easy settings by touch key operation
- 3 Graphic LCD screen
- 4 Pressure and system self-testing
- 5 Self-air pressure sensor calibration
- 6 Safety alarm system (cuff, pressure, battery, time check alarms)
- 7 Built-in battery
- 8 Stand-by mode
- 9 2 channels (4 ports for 2 double or single cuffs use)
- 10 RIVA feature
- 11 LOP sensor application
- 12 Optional 14 different cuff sizes (Single, IVRA, Conic)

Order Code: M067803

Product No: 8698864012745



C60

Specialized Neonatal Monitor

promedic
by **COMEN**

- 1 8.4" TFT touchscreen
- 2 Special design for neonatals
- 3 Up to 4 hours operating time with battery
- 4 Audible and visual alarms
- 5 At least 120 hours trend recording
- 6 Apnea alarm
- 7 Masimo SPO2 technology
- 8 Optional Apnea Self -saving monitoring
- 9 Optional Oxygen concentration monitoring function
- 10 Standard parameters: ECG, SPO2, RESP, TEMP, NIBP
- 11 Optional parameters: EtCO2, IBP, Thermal printer



Order Code: M066872
Product No: 8698864013698

NC3

Vital Signs Monitor

promedic
by **COMEN**

- 1 LED display
- 2 Suitable for adult, children and neonatal patients
- 3 Quick and accurate SPO2 and NIBP measurements
- 4 Optional TEMP module
- 5 Up to 12 hours operating time with battery



Order Code: M066870
Product No: 8698864013674

HECG

ECG (RHYTHM) HOLTER SYSTEM / Système Holter ECG (Rythme)
MULTIPLE LANGUAGE SELECTION / Sélection de plusieurs langues

- 1 Lightweight and compact design
- 2 Pacemaker analysis
- 3 Patient event marking
- 4 Use with 1.5 V AAA batteries
- 5 Advanced analysis by pc software
- 6 Turkish language
- 7 Optional: HL7 and DICOM connections



HECG-3

Order Code: M067406

Product No: 8698864013247

HECG-12

Order Code: M067401

Product No: 8698864013261

HOLTER SİSTEMİ YAZILIMI

Order Code: M067404

Product No: 8698864012660

HOSPICARE

Handheld Pulseoximeter

- 1 Measures blood oxygen saturation level and pulse rate
- 2 2.4" LCD display
- 3 Suitable for Adult, Pediatric & Neonatal
- 4 Audible and visual alarms
- 5 Low battery alarm
- 6 Patient data storage
- 7 Demo mode
- 8 Operating up to 30 hours with AA batteries



Order Code: M066606

Product No: 8698864014176

FAST & SIMPLE ECG BELTS

- 1 Intended for 12-lead Rest ECG
- 2 Compatible with all 12-lead ECG brands
- 3 Connect to the patient cable with snaps or crocodiles
- 4 Include 6 chest Ag/AgCl electrodes
- 5 With Handles or Strap
- 6 No maintenance
- 6 Electrodes and Plastic Parts may be replaced on-site



ECG BELT – With handles

Order Code: M061475

Product No: 8698864014183

ECG BELT – With Straps

Order Code: M061476

Product No: 8698864014190

LX

ELECTROSURGICAL UNIT – 200 WATT

- 1 6SENSE Tissue Feedback Technology
- 2 Self Control and Error Code Technology
- 3 Wide Range of Surgical Application
- 4 2 Bipolar, 6 Monopolar Modes
- 5 1 Monopolar, 1 adet Bipolar Output
- 6 1 User Adjustable Program
- 7 Patient Return Electrode Monitoring
- 8 Upgradability: Beam Coagulation / Ultrasonic Aspirator/ Smoke Evacuator
- 9 Audible and Visual Alarms

Order Code: M066851

Product No: 8698864013643



EASY +

ELECTROSURGICAL UNIT – 300 WATT

- 1 6SENSE Tissue Feedback Technology
- 2 Self Control and Error Code Technology
- 3 Wide Range of Surgical Application
- 4 2 Bipolar, 7 Monopolar Modes
- 5 Endoscopic (endocut) mode
- 6 1 Monopolar, 1 Bipolar Output
- 7 10 User Adjustable Programs
- 8 Patient Return Electrode Monitoring
- 9 Upgradability: Beam Coagulation / Ultrasonic Aspirator/ Smoke Evacuator
- 10 Audible and Visual Alarms

Order Code: M066852

Product No: 8698864013650



MX

ELECTROSURGICAL UNIT – 300 WATT

- 1 6SENSE Tissue Feedback Technology
- 2 Self Control and Error Code Technology
- 3 Wide Range of Surgical Application
- 4 4 Bipolar, 7 Monopolar Modes
- 5 Endoscopic (endocut) mode
- 6 2 Monopolar, 1 Bipolar Output
- 7 Simultaneous Coagulation
- 8 30 User Adjustable Programs
- 9 Patient Return Electrode Monitoring
- 10 Upgradability: Beam Coagulation / Ultrasonic Aspirator/ Smoke Evacuator
- 11 Audible and Visual Alarms



Order Code: M066850

Product No: 8698864013636

SIGMA +

ELECTROSURGICAL UNIT – 400 WATT

- 1 4,7" Colour LCD Display
- 2 6SENSE Tissue Feedback Technology
- 3 Self Control and Error Code Technology
- 4 Wide Range of Surgical Application
- 5 4 Bipolar, 8 Monopolar Modes
- 6 Endoscopic (endocut) mode
- 7 2 Monopolar, 1 Bipolar Output
- 8 Simultaneous Coagulation
- 9 30 User Adjustable Programs
- 10 Patient Return Electrode Monitoring
- 11 Upgradability: Beam Coagulation / Ultrasonic Aspirator/ Smoke Evacuator
- 12 Audible and Visual Alarms



Order Code: M066853

Product No: 8698864013667

STAR5000

Fetal & Maternal Monitor

- 1 12.1" Color TFT LCD touch screen
- 2 Advanced ARM system, convenient for maintenance and upgrade
- 3 Standard and Big Font parameter interfaces
- 4 High sensitivity waterproof probes
- 5 Built-in 112mm width thermal printer
- 6 Full data review up to 120 hours
- 7 Automatic fetal status scoring system
- 8 Wall-mountable
- 9 Wired/Wireless central monitoring system options up to 128 beds
- 10 Option of connectivity to external printer
- 11 3 levels alarm system
- 12 Standard configurations: FHR, TOCO, FM, AFM, Thermal printer
- 13 Optional configurations: FHR2, Fetal Stimulator, External Printer



Order Code: M066201

Product No: 8698864013629

STAR5000 E

Fetal & Maternal Monitor

- 1 5.7 inch TFT touchscreen
- 2 Built-in 112mm width thermal printer
- 3 Built-in Li-ion battery, up to 4 hours operating time
- 4 Manual / Automatic FM measurement
- 5 Automatic fetal status scoring system
- 6 FHR signal strength indicator
- 7 60 hours real-time waveform reply & storage
- 8 Timing monitoring function to avoiding over-time monitoring
- 9 Optional connectivity to central monitoring system



Order Code: M066200

Product No: 8698864013612

CLARMAC/CLARMIL

STANDARD LARYNGOSCOPE SET
MACINTOSH/MILLER

- 1 Within plastic case
- 2 Autoclavable at 134°C
- 3 Curved and straight blade types
- 4 Resistant to disinfectant solution
- 5 The set is composed of:
 - o 3 pieces of blade (Sizes: 00, 0, 1, 2, 3, 4)
 - o 1 piece of handle with special lock system
 - o 1 piece of spare bulb
 - o Optional in sets of 3, 4, 5 blades



CLARMAC

Order Code:

Product No: 8698864008595

CLARMIL

Order Code:

Product No: 8698864008625

FLARMAC/FLARMIL

FIBER OPTIC LARYNGOSCOPE SET
MACINTOSH/MILLER

- 1 Within plastic case
- 2 Autoclavable at 134°C
- 3 Curved and straight blade types
- 4 Resistant to disinfectant solution
- 5 The set is composed of:
 - o 3 pieces of blade (Sizes: 00, 0, 1, 2, 3, 4)
 - o 1 piece of handle with special lock system
 - o 1 piece of spare bulb
 - o Optional in sets of 3, 4, 5 blades



FLARMAC

Order Code:

Product No: 8698864008649

FLARMIL

Order Code:

Product No: 8698864008663

DIASET

DIAGNOSTIC SET
STANDARD/HALOGEN

- 1 Within zipper case
- 2 Autoclavable at 134°C
- 3 Adjustable light intensity
- 4 The set is composed of:
 - o 1 piece of otoscope
 - o 1 piece of ophthalmoscope
 - o 2 pieces of larynx mirror (sizes of 3 and 4)
 - o 1 piece of metal tongue depressor
 - o 1 piece of adapter for disposable tongue depressor
 - o 3 pieces of ear speculum (sizes of 2,5/3,5/4,5 mm)
 - o 1 piece of light holder (for larynx mirrors)
 - o 2 pieces of auxiliary bulb
 - o 1 piece of nose speculum (reusable)
 - o 1 piece of handle with screw lock system



Order Code: M040410

Product No: 8698864008687

OTOSET

OTOSCOPE SET / Ensemble d'Otoscope
STANDARD/HALOGEN

- 1 Within zipper pouch
- 2 Adjustable light intensity
- 3 The set is composed of:
- 4
 - o 1 piece of otoscope head with screw lock system
 - o 3 pieces of reusable ear speculum (2,5/3,5/4,5 mm)
 - o 1 piece of handle with screw lock system



Order Code: M040411
Product No: 8698864008700

OFTASET

OPHTALMASCOPE SET / Ensemble d'Ophthalmoscope
STANDARD/HALOGEN

- 1 Within plastic case
- 2 Autoclavable at 134°C
- 3 24 diopter setting between -25 and +40
- 4 The set is composed of:
 - 1 piece of ophtalmoscope head with screwed lock system
 - 1 piece of handle with screw lock system



Order Code: M040412
Product No: 8698864008700

OTO-OFTASET

OTO/OPHTALMASCOPE SET
Ensemble de Otoscope + Ophthalmoscope
STANDARD/HALOGEN

- 1 Within plastic case
- 2 Adjustable light intensity
- 3 Vacuum type bulb
- 4 24 diopter setting at a range between -25 and +40
- 5 The set is composed of:
 - o 3 pieces of reusable ear speculum (2.5/3.5/4.5 mm)
 - o 1 piece of otoscope head
 - o 1 piece of ophthalmoscope head
 - o 1 piece of handle with screw lock system



Order Code: M040413
Product No: 8698864008724



firstMED



Index



PATIENT MONITORS.....	50
ECG UNITS	51 - 52
ECG HOLLERS.....	53
BLOOD PRESSURE HOLLER.....	53
SUCTION UNITS.....	54
OXYGEN REGULATOR.....	54
INCUBATORS.....	55
SPIROMETRY SYSTEMS.....	55

pM-6000

Patient Monitor - Compact Design

- 1 8.4" LCD color screen
- 2 2 hours of battery life
- 3 Simultaneous 8 waveforms displaying
- 4 Visual and audible alarms
- 5 72 hours of trend saving
- 6 Standard parameters: ECG, SPO2, TEMP, NIBP, HR, RESP
- 7 Optional parameters: 2 Channels IBP, CO2

Order Code: M068101

Product No:8698864008380



pM-6800

Patient Monitor - Compact Design

- 1 8" LCD color screen
- 2 2 hours of battery life
- 3 Simultaneous 7 waveforms displaying
- 4 Visual and audible alarms
- 5 72 hours of trend saving
- 6 Apnea alarm
- 7 Standard parameters: ECG, SPO2, TEMP, NIBP, HR, RESP

Order Code: M068102

Product No:8698864008397



pM-8000

Patient Monitor - Compact Design

- 1 12.1" LCD color screen
- 2 2 hours of battery life
- 3 Simultaneous 8 waveforms displaying
- 4 Visual and audible alarms
- 5 72 hours of trend saving
- 6 Apnea alarm
- 7 Standard parameters: ECG, SPO2, TEMP, NIBP, HR, RESP
- 8 Optional parameters: 2xIBP, ETCO2

Order Code: M068103

Product No:8698864008373



ECG-100

ECG Device
Single Channel / Tabletop

- 1 Lightweight and compact design
- 2 Print modes: Automatic and manual
- 3 Recording system: Thermal printer
- 4 AC and EMG filters
- 5 Internal battery

Order Code: M066922

Product No:8698864008915



ECG-300

ECG Device
3 Channels / Tabletop

- 1 Lightweight and compact design
- 2 Print modes: Automatic, manual and rhythm
- 3 Recording system: Thermal printer
- 4 AC, EMG and DFT filters
- 5 Internal battery
- 6 Internal memory: 1.000 patients
- 7 Lead placement diagram
- 8 Display: 320x240 color LCD

Order Code: M066923

Product No:8698864008939



ECG-600

ECG Device

6 Channels/Tabletop/Touch Screen

- 1 Lightweight and compact design
- 2 Print modes: Automatic, manual and rhythm
- 3 Recording system: Thermal printer
- 4 AC, EMG and DFT filters
- 5 Internal memory: 1.000 patients
- 6 Lead placement diagram
- 7 Display: 800x480 color LCD touch screen
- 8 Interpretation function



Order Code: M066925

Product No:8698864008946

ECG-1200

ECG Device

12 Channels/Tabletop/Touch Screen

- 1 Lightweight and compact design
- 2 Print modes: Automatic, manual and rhythm
- 3 Recording system: Thermal printer
- 4 AC, EMG and DFT filters
- 5 Internal battery
- 6 Internal memory: 1.000 patients
- 7 Lead placement diagram
- 8 Display: 800x600 color LCD touch screen



Order Code: M066926

Product No:8698864008953

EH-12

ECG Holter System
12 Channels

- 1 Lightweight and compact design
- 2 12 channels ECG recording
- 3 1.2 GB built-in memory
- 4 Analysis channel selection
- 5 Patient event marking
- 6 Template, ST, HRV, QT, Atrial Fibrillation, HRT and T-Wave analyzes
- 7 Advanced reporting

Order Code: M067121

Product No: 8698864008885



ABPM-50

Ambulatory Blood Pressure Holter System
Large and Color Screen

- 1 Compact and lightweight design
- 2 Oscillometric measurement
- 3 Windows based PC software
- 4 Programmable software for day and night
- 5 Advanced reporting
- 6 Audible and visual alarms

Order Code: M067102

Product No: 8698864012608



SP-9

Suction Unit

- 1 Suitable for emergency service, ambulance, first aid and hospital usage
- 2 No oil and maintenance free motor
1 liter plastic jar
- 3 Maximum negative pressure: 0.075 MPa
- 4 Power supply: AC 220V or DC 12V
- 5 Vacuum rate: 15 lt/min and 20 lt/min

Order Code: M067308

Product No:8698864013063



OXYFLOW-01

Oxygen Regulator
Medical Use

- 1 Oxygen flow range: 1 - 15L/min
- 2 Pressure reducing range: 2.0 - 3.0 kg
- 3 Relief valve pressure range: (3.5 kg \pm 0.5 kg)

Order Code: M068250

Product No:8698864011274



INC-100

Incubator
Microprocessor Controller

- 1 Control mode: Air mode
- 2 Two level humidity adjustment
- 3 6 operating windows
- 4 RS232 connector
- 5 Optional parts: Up phototherapy unit and down phototherapy unit

Order Code: M068220
Product No: 8698864013087



SP-10

Spirometry System
Portable


- 1 1.8" TFT color screen
- 2 Test items: FVC
- 3 Internal battery

Order Code: M066910
Product No: 8698864009325





PLASMED



Index



SURGICAL BLADES	57
RESUSCITATORS	57
SKIN STAPLER & REMOVER.....	58
ELECTROSURGICAL PADS.....	58
INFUSION SETS	58
I.V. CANNULA	59
ECG ELECTRODE	59
LATEX EXAMINATION GLOVES	60
LATEX SURGICAL GLOVES.....	60

SBL

SURGICAL BLADE / DISPOSABLE / STERILE



- 1 Carbon steel
- 2 Sizes:
10, 11, 12, 15, 18, 20, 21, 22, 23, 24
- 3 Individual sterile package
- 4 Packing: 100 pcs/box, 5.000 pcs/carton

SBL-10
Order Code: M035011
Product No: 6945630106415

SBL-20
Order Code: M035016
Product No: 6945630106460

SBL-11
Order Code: M035012
Product No: 6945630106422

SBL-21
Order Code: M035017
Product No: 6945630106477

SBL-12
Order Code: M035013
Product No: 6945630106439

SBL-22
Order Code: M035018
Product No: 6945630106484

SBL-15
Order Code: M035014
Product No: 6945630106446

SBL-23
Order Code: M035019
Product No: 6945630106491

SBL-18
Order Code: M035015
Product No: 6945630106453

SBL-24
Order Code: M035020
Product No: 6945630106507

AMBUSET-DPVC RESUSCITATOR KIT
DISPOSABLE

- 1 Pressure relieving valve
- 2 Inflatable aircushion mask with syringe
- 3 Oxygen balloon volumes: Adult (1700ml), pediatric (600ml) and infant (280ml)
- 4 Oxygen reservoir volumes: Adult (2000 ml), pediatric and infant (1600 ml)
- 5 1 piece PVC mask, 3 pieces airway and 1 piece nasal oxygen cannula in kit
- 6 Packing: 24 pcs/carton



D280
Order Code: M034523
Product No: 8698864013193

D600
Order Code: M034524
Product No: 8698864013209

D1700
Order Code: M034525
Product No: 8698864013216

AMBUSET-RSILICONE RESUSCITATOR KIT
REUSABLE

- 1 Autoclavable (134°C)
- 2 Pressure relieving valve
- 3 Oxygen balloon volumes: Adult (1700ml), pediatric (600ml) and infant (280ml)
- 4 Oxygen reservoir volumes: Adult (2000 ml), pediatric and infant (1600 ml)
- 5 Mask sizes: 0, 1, 2, 3, 4, 5
- 6 2 pieces silicone mask, 3 pieces airway, 1 piece nasal oxygen cannula and 1 piece mouth opener in kit
- 7 Hard plastic carrying case in order to carry all the equipments
- 8 Packing: 12 pcs/carton



R280
Order Code: M034501
Product No: 8698864013162

R600
Order Code: M034521
Product No: 8698864013179

R1700
Order Code: M034522
Product No: 8698864013186

F-35W

STAPLER REMOVER / REUSABLE / STERILE

- 1 Can be sterilized by ethylene oxide (ETO)
- 2 Individual sterile blister package
- 3 Packing: 200 pcs/carton



Order Code: M035031
Product No: 5060186750125

C-WR

SKIN STAPLER / DISPOSABLE / STERILE

- 1 Stainless steel
- 2 35 pcs wide (w) hooks
- 3 Individual sterile blister package
- 4 Packing: 5 pcs/box, 100 pcs/carton



Order Code: M035030
Product No: 5060186751450

GP03 / GP04

ELECTROSURGICAL PAD
DISPOSABLE

- 1 Monopolar or bipolar
- 2 Hydrogel conductive adhesive
- 3 Compatible with all electrosurgical units
- 4 Reusable cables (Optional)
- 5 Packing: 50 pcs/box, 400 pcs/carton



GP03
Order Code: M034302
Product No: 8698864001404

GP04
Order Code: M034304
Product No: 8698864001428

INFSET-S

INFUSION SET (VENTED) DISPOSABLE / STERILE

- 1 For infusion bottles and bags
- 2 Vented spike
- 3 Air-inlet filter with snap on cap
- 4 20 drops= 1 ml
- 5 15 µm filter
- 6 Roller clamp
- 7 injection site
- 8 Luer slip connector
- 9 21 G needle
- 10 Tube length: 150 cm
- 11 Individual sterile package
- 12 Packing: 500 pcs/carton



Order Code: M034231
Product No: 8698864009912

PLASCATH

I.V. CANNULA

- 1 Japanese Stainless Steel Cannula
- 2 Extra Fine Needle Point for Smooth Cannulation
- 3 Thin Walled FEP/PTFE Catheter for Lower Penetration Force
- 4 Customized Tipping Technology for Lower Tip-Wall Thickness
- 5 Minimum Clearance between Catheter and Needle for Zero Peel-Back Issues during use
- 6 Luer Lock Designed to Prevent Leakage
- 7 Non-pyrogenic
- 8 No: 14G, 16G, 17G, 18G, 20G, 22G, 24G With Wings & Injection Port Radio-opaque Catheter for Adult
24G, 26G With Small Wings & Without Injection Port Radio-opaque Catheter for Neonatal

No:14G Adult
Order Code: M031825
Product No: 8698864014671

No:16G Adult
Order Code: M031826
Product No: 8698864014688

No:17G Adult
Order Code: M031827
Product No: 8698864014695

No:18G Adult
Order Code: M031828
Product No: 8698864014701


No:20G Adult
Order Code: M031829
Product No: 8698864014718

No:22G Adult
Order Code: M031830
Product No: 8698864014725

No:24G Adult
Order Code: M031831
Product No: 8698864014732

No:24G Neonatal
Order Code: M031841
Product No: 8698864014749

No:26G Neonatal
Order Code: M031842
Product No: 8698864014756




ECG ELECTRODE

Available Types

Adult: Round, Diameter 50 mm, Foam Holder, Solid Gel, 50 pcs / pouch.

Adult: Round, Diameter 50 mm, Non-Woven Bracket, Liquid Gel, 50 pcs / pouch.

Child: 32 x 38 mm Foam Holder, Solid Gel, 50 pcs / pouch.



ELEC-S55
Order Code: M035920
Product No: 8698864013407

ELEC-S30
Order Code: M035921
Product No: 8698864013414

ELEC-W55
Order Code: M035922
Product No: 8698864013421

LATEX EXAMINATION GLOVES

- 1 Ambidextrous
- 2 Single use only
- 3 Powdered / Powder-free
- 4 Sizes: X-Small, Small, Medium, Large, X-Large



POWDERED / X-SMALL

Product No:
8698864013834

POWDERED / SMALL

Product No:
8698864013841

POWDERED / MEDIUM

Product No:
8698864013858

POWDERED / LARGE

Product No:
8698864013865

POWDERED / X-LARGE

Product No:
8698864013872

POWDERED-FREE / X-SMALL

Product No:
8698864013889

POWDERED-FREE / SMALL

Product No:
8698864013896

POWDERED-FREE / MEDIUM

Product No:
8698864013902

POWDERED-FREE / LARGE

Product No:
8698864013919

POWDERED-FREE / X-LARGE

Product No:
8698864013926

LATEX SURGICAL GLOVES

- 1 Natural rubber latex
- 2 Single use only
- 3 Powdered / Powder-free
- 4 Sizes : 6 - 6,5 - 7 - 7,5 - 8 - 8,5



POWDERED / No: 6

Order Code: M010160
Product No:8698864014823

POWDERED / No: 6,5

Order Code: M010161
Product No:8698864014770

POWDERED / No: 7

Order Code: M010162
Product No:8698864014787

POWDERED / No: 7,5

Order Code: M010163
Product No:8698864014794

POWDERED / No: 8

Order Code: M010164
Product No:8698864014800

POWDERED / No: 8,5

Order Code: M010165
Product No:8698864014817

POWDERED-FREE / No: 6

Order Code: M010166
Product No:8698864014824

POWDERED-FREE / No: 6,5

Order Code: M010167
Product No:8698864014831

POWDERED-FREE / No: 7

Order Code: M010168
Product No:8698864014838

POWDERED-FREE / No: 7,5

Order Code: M010169
Product No:8698864014845

POWDERED-FREE / No: 8

Order Code: M010170
Product No:8698864014862

POWDERED-FREE / No: 8,5


Order Code: M010177
Product No:8698864014879



DISTRIBUTION PRODUCTS



Index




MEDCAPTAIN.....	62 - 63
CHEST.....	63
KAI.....	64
NULIFE.....	64
CA-MI.....	65
F.BOSCH.....	65

MP-60

INFUSION PUMP - MODULAR DESIGN

MEDCAPTAIN

- 1 Touchscreen
- 2 Lightweight and small
- 3 Modular desing
- 4 Plug&play
- 5 Stackable
- 6 Automatic door
- 7 Anti free flow mechanism
- 8 Automatic and manual bolus
- 9 Adjustable KVO
- 10 Night mode
- 11 Patient data
- 12 Drug library
- 13 IV standına sabitleme imkânı
- 14 Fixing to IV pole
- 15 Audible and visual alarms
- 16 Internal battery



Order Code: M066830
Product No:8698864012639

SYS-6010

INFUSION PUMP – PORTABLE

MEDCAPTAIN

- 1 User-friendly
- 2 Touchscreen
- 3 Anti free flow mechanism
- 4 High safety and stability
- 5 High accuracy
- 6 Optional drop sensor




Order Code: M066831
Product No:8698864012646

MP-80

INFUSION WORKSTATION

MEDCAPTAIN

- 1 Modular desing
- 2 Stackable
- 3 Plug&Play
- 4 Synchronize alarm with pumps
- 5 Wi-Fi Connectivity to central monitoring system
- 6 Compatible with MP-60 infusion pump and MP-30 syringe pump




Order Code: M066834
Product No:8698864012653

MP-30

SYRINGE PUMP - MODULAR DESIGN

MEDCAPTAIN

- 1 Touchscreen
- 2 Lightweight and small
- 3 Modular desing
- 4 Plug&play
- 5 Stackable
- 6 Automatic and manual bolus
- 7 Adjustable KVO
- 8 Night mode
- 9 Patient data
- 10 Drug library
- 11 Fixing to IV pole
- 12 Audible and visual alarms
- 13 Internal battery



Order Code: M066832
Product No:8698864012615

SYS-3011

SYRINGE PUMP
PORTABLE

MEDCAPTAIN

- 1 User-friendly
- 2 High safety and stability
- 3 High accuracy
- 4 Patented drive system
- 5 Special technology of motor drive



Order Code: M066833
Product No:8698864012622

HI-101

SPIROMETRY SYSTEM
TABLE TOP

CHEST

- 1 For testing asthma, COPD and other respiratory failure patients
- 2 Compact for mobility
- 3 Large LCD screen
- 4 Test items: SVC, FVC, MVV, BD
- 5 Proven ATS accuracy
- 6 Data storage: Approx. 100 tests
- 7 Computer interface: RS-232C
- 8 Ideal for clinics
- 9 Just pres the buttons and testing is completed in 3 minutes
- 10 Rechargeable lead acid battery
- 11 Thermal printer
- 12 Data management software for Windows
- 13 Optional items: Calibration syringe (3L or 1 L), spiro filtre



Order Code: M063901
Product No:8698864013353

BIOPSY PUNCH

STERILE

- 1 Disposable
- 2 Dimensions: 2 mm, 3 mm, 4 mm, 5 mm, 6 mm
- 3 Packing: 20 pcs/box; 50 boxes/carton



2 mm

Order Code: M032307

Product No:4560146922605

3 mm

Order Code: M032301

Product No:4560146922599

4 mm

Order Code: M032302

Product No:4560146922582

5 mm

Order Code: M032303

Product No:4560146922612

6 mm

Order Code: M032304

Product No:4560146922575

Male External (Condom) Catheter

DISPOSABLE / STERILE

- 1 Soft, pure latex
- 2 Compatibility to all standard urine bags
- 3 With soft adhesive tapes
- 4 Sizes: Medium (25 mm), large (30 mm) and extra large (35 mm)
- 5 Individual sterile blister package
- 6 Packing: 50 pcs/box, 600 pcs/carton

CATHETER-L

Order Code: M030702

Product No:8908000123670

CATHETER-M

Order Code: M030701

Product No:8908000123663

CATHETER-XL

Order Code: M030703

Product No:8908000123687



NEW ASKIR 230/12V BR

BATTERY-OPERATED SUCTION UNIT

NEW ASKIR 230/12V BR is an electric suction unit for the aspiration of body liquids, oral, nasal and tracheal aspiration in adults or children. Suitable for emergency applications and for post-operative therapy, both in professional environment and home-care.

Three different options for operation: AC/DC adapter, rechargeable battery and 12V make it a versatile suction unit approved for NON-STOP operation without overheating. Provided with visual and acoustic alarm indicating low battery level.

The newly-designed ergonomic handle ensures a comfortable and secure grip.

Order Code: M066933
Product No:8054619100370



KONSTANTE II

Aneroid Sphygmomanometer
Professional



F. BOSCH
MEDIZINTECHNIK
GERMANY

- 1 Diameter scale: 50 mm (0-300 mmHg)
- 2 Double tubing system
- 3 Soft carrying case
- 4 Hook cuff (Velcro cuff-Optional)
- 5 Adult and children size cuffs



ADULT

Order Code: M020301
Product No:4260157891012

CHILD

Order Code: M020307
Product No:4260157891029

Various TRIMPEKS products which do not take place in the Product Catalogue

Hospital Furniture

Hospitalisation Bed
Baby Cot
Mayo Table
Nightstand
Delivery Table
Operating Table
Examination Table
Blood Sampling Chair
Dialysis Chair
Ward Screen
Couch Steps
Stretcher
Trolleys / Carts
Mobile Infusion Stand (Serum holder)
Adjustable Stool
Film Viewer
Surgical Sink
Hospital Waiting Bench

Miscellaneous

ENT Consultation Unit
Dialysis Machine
Medical Refrigerator
ICU Ventilator
Emergency & Transport Ventilator
Anesthesia Machine
Defibrillator
AED (Automated External Defibrillator)
Endoscopy Column
Laparoscopy Set
Feces Analyzer
Dental Chair

Radiology Products

Mobile X-Ray
Ultrasound
Dental Panoramic X-Ray
CT Scanner

Lamps

Examination Lamp
Surgical Lamp

Digital and Mechanical Scales

Digital Baby Scale
Mechanical Baby Scale with sliding weights
Mechanical Adult Scale with sliding weights and Height Rod

Consumables

IV Cannula
Infusion Set
Disposable Syringes
Angiographic Syringes
Plasters, Bandages, Compresses, Dressings, Gauze Rolls, Cotton, etc.

RUTF (Ready to Use Therapeutic Food)

Produits divers de TRIMPEKS qui ne figurent pas dans le Catalogue de Produits

Mobilier Hospitalier

Lit d'hospitalisation

Berceau

Table de Mayo

Table de chevet

Table d'accouchement

Table d'operation

Chaise / Lit / Table / Divan d'Examen

Fauteuil de prelevement

Fauteuil de dialyse

Paravent

Escabeau

Brancard

Chariots

Porte-sérum (Potence)

Tabouret réglables

Négatoscope

Lavabo chirurgical

Banc d'attente

Radiologie

Radiographie Mobile

Ecographe

Radio Panoramique

Scanner CT

Lampes

Lampe d'examen

Lampe scialytique

Balances numériques et mécaniques

Pèse-bébé numérique

Pèse bébé mécanique à poids

coulissant

Pèse personne mécanique à poids

coulissant

Divers

Unité de consultation ORL

Generateur de dialyse

Réfrigérateur médicalisé

Ventilateur de soins intensifs, Ventilateur d'urgence et de transport

Appareil d'anesthésie

Défibrillateur

Appareil AED

Colonne d'endoscopie

Ensemble de coelioscopie

Analyseur fécal

Fauteuil dentaire

Consommables

Cathéter intraveineuse

Set de perfusion intraveineuse

Seringues jetables

Seringues angiographiques

Sparadraps, Rouloaux de Gaze,

Bandages, Compresses,

Pansements, Cotton, etc.,

ATPE (Aliments thérapeutiques prêts à l'emploi)

TRIMPEKS

Your
“One-Stop Shop”
for Healthcare
Requirements for

- ✓ **Detection**
- ✓ **Risk Assessment**
- ✓ **Diagnosis**
- ✓ **Therapy**



Contact Us

 Sultan Selim Mah. Yunus Emre Cd. No:1/11 34415 Kağıthane İstanbul Türkiye

 www.trimpeks.com

 export@trimpeks.com

 +90 212 319 50 00